

---

**Product Data Sheet**

---

Product Name: Rigin  
Cat. No.: GA23485

**Chemical Properties**

Cas. No. 77727-17-4

Formula  $C_{18}H_{32}N_8O_6$  M.Wt 456.5

Solubility Soluble in DMSO Storage Store at  $-20^{\circ}C$ , protect from light

General tips For obtaining a higher solubility, please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

**Background**

Besides tuftsin, rigin is another phagocytosis-stimulating tetrapeptide originally isolated from human IgG.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA