
Product Data Sheet

Product Name: Sar-Pro-Arg-pNA

Cat. No.: GA23495

Chemical Properties

Cas. No. 77695-30-8

Formula $C_{20}H_{30}N_8O_5$ M.Wt 462.51

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**Chromogenic substrate for the assay of human α -thrombin. $K_m = 75 \mu M$.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA