

---

**Product Data Sheet**

---

Product Name: SASRIN™ resin (100-200 mesh, 0.7-1.0 mmol/g)

Cat. No.: GA23497

**Chemical Properties**

Cas. No. 124760-64-1

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

**Background**

For further details see SASRIN resin (200-400 mesh, 0.80-1.20 mmol/g).

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA