
Product Data Sheet

Product Name: SPARC (119-122) (mouse)

Cat. No.: GA23510

Chemical Properties

Cas. No. 155149-79-4

Formula $C_{20}H_{36}N_8O_5$ M.Wt 468.56

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

SPARC is a transiently expressed extracellular matrix-binding protein that alters cell shape and regulates endothelial cell proliferation in vitro. KGHK, a copper-binding peptide fragment, has been shown to be a potent stimulator of angiogenesis. Speci

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA