
Product Data Sheet

Product Name: Stresscopin (human)

Cat. No.: GA23515

Chemical Properties

Cas. No. 352020-03-2

Formula $C_{195}H_{326}N_{56}O_{53}S_2$

M.Wt 4367.21

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Stresscopin is a specific ligand for the corticotropin-releasing hormone receptor type 2 (CRHR2). Transcripts are expressed in brain and most tissues analyzed. Intraperitoneal injections of stresscopin suppresses heat-induced edema formation in anesthe

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA