

---

**Product Data Sheet**

---

Product Name: Suc-Ala-Ala-Pro-Ala-pNA

Cat. No.: GA23527

**Chemical Properties**

Cas. No. 72682-69-0

Formula  $C_{24}H_{32}N_6O_9$  M.Wt 548.55Solubility Soluble in DMSO Storage Store at  $-20^{\circ}C$ 

General tips For obtaining a higher solubility, please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Suc-AAPA-pNA, an excellent substrate for porcine pancreatic elastase. Substrate for  $\alpha$ -lytic protease.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA