
Product Data Sheet

Product Name: Suc-Ala-Ala-Pro-Arg-pNA

Cat. No.: GA23528

Chemical Properties

Cas. No. 131068-47-8

Formula C₂₇H₃₉N₉O₉ M.Wt 633.66

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Suc-AAPR-pNA, a substrate for wild-type and K188D/D189K trypsins (K_m = 32.8 mM respectively 50.8 mM).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA