
Product Data Sheet

Product Name: Suc-Ala-Leu-Pro-Phe-pNA

Cat. No.: GA23547

Chemical Properties

Cas. No. 128802-78-8

Formula C₃₃H₄₂N₆O₉ M.Wt 666.73

Solubility DMSO : ≥ 100 mg/mL (149.99 mM) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Substrate for FK-506 binding proteins (FKBPs, also called macrophilins) and cyclophilins, which belong to the group of peptidyl prolyl cis-trans isomerases (PPIases). Suc-ALPF-pNA, has been used for an uncoupled protease-free assay of PPIase activity.</p></div>

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
