

---

**Product Data Sheet**

---

Product Name: Tau Peptide (244-274) (Repeat 1 Domain)

Cat. No.: GA23577

**Chemical Properties**

Cas. No. 330456-24-1

Formula  $C_{141}H_{238}N_{42}O_{44}S$ 

M.Wt 3257.63

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

During heparin-induced self-aggregation of the four repeat domain peptides (R1-R4) excised from tau, R1, which resembles Tau Peptide (244-274), was the most resistive to filament formation. Its conformation remained unchanged. The peptide maintained th

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA