
Product Data Sheet

Product Name: Tau Peptide (294-305) (human)

Cat. No.: GA23585

Chemical Properties

Cas. No. 1428135-29-8

Formula $C_{51}H_{85}N_{17}O_{17}$ M.Wt 1208.32

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Tau 294-305, KDNKHVPGGGS, conjugated to KLH via an N-terminal cysteine was tested as anti-tau vaccine. The tau fragment might prevent formation of the beta/beta-helix structure of protofilaments involving residues 306-378.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA