
Product Data Sheet

Product Name: Tau Peptide (323-335)

Cat. No.: GA23593

Chemical Properties

Cas. No. 2022956-59-6

Formula $C_{60}H_{103}N_{21}O_{20}$ M.Wt 1438.61

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The tau epitope was used for raising the anti-tau antibody BR135 directed against the microtubule-binding repeat region.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA