
Product Data Sheet

Product Name: Tau Peptide (379-408)

Cat. No.: GA23595

Chemical Properties

Cas. No. 1246357-53-8

Formula C₁₄₀H₂₂₈N₄₄O₄₆

M.Wt

3263.62

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The non-phosphorylated analog of the tau fragment 379-408 phosphorylated at Ser³ and Ser could partially block the binding of antibodies against the 396/404 tau region.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA