
Product Data Sheet

Product Name: Tau Peptide (45-73) (Exon 2/Insert 1 Domain)

Cat. No.: GA23596

Chemical Properties

Cas. No. 2022956-54-1

Formula $C_{120}H_{189}N_{31}O_{57}$

M.Wt 2978

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Tau Peptide (45-73) (Exon 2/Insert 1 Domain) is a sequence of the so-called "projection domain" of tau. Phosphorylated residues 68, 69 and 71 of this sequence have been identified on the transformed PHF-tau (paired helical filament-tau).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA