
Product Data Sheet

Product Name: Thrombin Receptor Antagonist

Cat. No.: GA23607

Chemical Properties

Cas. No. 207553-92-2

Formula C₅₈H₉₁N₁₇O₂₀S₂

M.Wt

1410.6

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The peptide T1, MSRPACPNDKYE, is based on the sequence of the thrombin receptor domain that functions as a tethered ligand for the receptor itself. T1 has been selected from a phage peptide library as an inhibitor of thrombin triggered platelet aggrega

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA