
Product Data Sheet

Product Name: Toolbox Fmoc-Amino Acid Derivatives

Cat. No.: GA23618

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The Toolbox contains the 20 standard Fmoc-amino acid L-derivatives at 5 g per each vial. These are Fmoc-Ala-OH . HO (No. 4003167), Fmoc-Arg(Pbf)-OH (No. 4089151), Fmoc-Asn(Trt)-OH (No. 4014948), Fmoc-Asp(OtBu)-OH (No. 4093142), Fmoc-Cys(Trt)-OH (No. 40

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA