

---

**Product Data Sheet**

---

Product Name: Tos-Gly-Pro-Lys-pNA

Cat. No.: GA23625

**Chemical Properties**

Cas. No. 73392-19-5

Formula  $C_{26}H_{34}N_6O_7S$  M.Wt 574.66

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Chromogenic serine protease substrate, e.g. for tryptase. Tos-GPK-pNA has been used for assaying bovine duodenase.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA