

---

## Product Data Sheet

---

Product Name: TRAF6 Peptide

Cat. No.: GA23635

### Chemical Properties

Cas. No. 852805-92-6

Formula  $C_{145}H_{238}N_{34}O_{44}$

M.Wt

3161.64

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

The cell-permeable peptide inhibits the interaction of TRAF6, an E3 ubiquitin ligase, with p62 (sequestosome 1). P62 has been shown to be a polyubiquitin shuttling factor which binds to K63-polyubiquitinated proteins. TRAF6-p62 interaction modulators m

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA