
Product Data Sheet

Product Name: Tricyclic amide linker resin (DL-form) (200-400 mesh, 0.5-0.8 mmol/g)
Cat. No.: GA23642

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Ramage resin is highly suitable for the preparation of peptide amides by the Fmoc-strategy and is recommended for C-terminal Phe, Tyr and Ile. Cleavage has been performed with 50 % TFA in CHCl₃ or 95% aqueous TFA. Scavengers may be required.
Tricyclic a

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA