
Product Data Sheet

Product Name: Tyr-(D-Dab⁴,Arg⁵,D-Trp⁸)-cyclo-Somatostatin-14 (4-11)

Cat. No.: GA23649

Chemical Properties

Cas. No. 496849-46-8

Formula C₆₇H₈₅N₁₅O₁₁ M.Wt 1276.51

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

High affinity ligand for all five somatostatin receptors (sst1-sst5). In AtT-20 mouse anterior pituitary tumor cells which mainly express sst2 and sst5 receptors the binding affinity of this somatostatin analog was found to be higher than that of somat

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA