

---

**Product Data Sheet**

---

Product Name: Tyr-LL-37

Cat. No.: GA23653

**Chemical Properties**

Cas. No. 2022972-74-1

Formula C<sub>214</sub>H<sub>349</sub>N<sub>61</sub>O<sub>55</sub>

M.Wt

4656.5

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

LL-37 analog for radioiodination.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA