

---

## Product Data Sheet

---

Product Name: Tyrosinase (192-200) (human, mouse)

Cat. No.: GA23654

### Chemical Properties

Cas. No. 170294-35-6

Formula  $C_{54}H_{77}N_{13}O_{17}$  M.Wt 1180.28

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

SEIWRDIDF, antigenic peptide presented by HLA-B44 and recognized on a human melanoma by autologous cytolytic T lymphocytes.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA