
Product Data Sheet

Product Name: Urocortin II (mouse)

Cat. No.: GA23659

Chemical Properties

Cas. No. 330648-32-3

Formula $C_{187}H_{320}N_{56}O_{50}$

M.Wt

4152.95

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Urocortin II (mouse) exhibited considerably higher affinity for both splice variants of the type 2 CRF receptor (CRF-R2) than for CRF-R1. The potencies of Ucn II (mouse) in activating CRF-R2 α and CRF-R2 β were nearly equivalent to that of uroc

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA