

---

**Product Data Sheet**

---

Product Name: Xanthenyl linker resin (200-400 mesh, 0.4-0.7 mmol/g)

Cat. No.: GA23669

**Chemical Properties**

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

**Background**

Sieber amide resin has been developed for the Fmoc-SPPS of fully t-butyl protected peptide amides. Cleavage can be achieved with 1% TFA in methylene chloride. Scavengers may be required. As N-methylation increases the acid-lability of the peptide bond,

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA