
Product Data Sheet

Product Name: Z-Ala-Pro-pNA

Cat. No.: GA23698

Chemical Properties

Cas. No. 66382-56-7

Formula $C_{22}H_{24}N_4O_6$ M.Wt 440.46

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The prolyl oligopeptidase (POP) substrate Z-AP-pNA was used in combination with Z-GP-pNA and Z-AA-pNA for characterizing the prolyl oligopeptide of of the hyperthermophile *Pyrococcus furiosus*.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA