

---

**Product Data Sheet**

---

Product Name: Z-Asp-Gln-Met-Asp-AFC

Cat. No.: GA23727

**Chemical Properties**

Cas. No. 1926163-45-2

Formula C<sub>36</sub>H<sub>39</sub>F<sub>3</sub>N<sub>6</sub>O<sub>13</sub>S

M.Wt 852.8

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Ac-DQMD-AFC, fluorogenic substrate for caspase-3. It contains the DQMD peptide sequence of the baculovirus caspase inhibitor p35.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---