
Product Data Sheet

Product Name: Z-Gly-Ala-Pro- β NA

Cat. No.: GA23793

Chemical Properties

Cas. No. 202001-62-5

Formula $C_{28}H_{30}N_4O_5$ M.Wt 502.57Solubility Soluble in DMSO Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Z-GAP- β NA was used in a staining assay for detecting thermostable mutants of prolyl endopeptidase expressed in E.coli.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA