
Product Data Sheet

Product Name: Z-Phe-Arg-4MβNA . HCl

Cat. No.: GA23893

Chemical Properties

Cas. No. 100900-17-2

Formula C₃₄H₃₈N₆O₅ · HCl

M.Wt 647.17

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Substrate for cathepsin B, cathepsin L, and cathepsin S. Inactivation of cathepsin B by 4M urea led to a selective assay of cathepsin L activity employing Z-FR-MNA in goat brain homogenates.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA