
Product Data Sheet

Product Name: α -Synuclein Binding Peptide

Cat. No.: GA24001

Chemical Properties

Cas. No.

Formula $C_{45}H_{80}N_{14}O_{14}$ M.Wt 1041.2

Solubility Soluble in DMSO Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Cheruvvara and coworkers conducted a single step selection followed by two rounds of competition selection protein-fragment complementation assay (PCA) in order to derive α -synuclein (α -syn) binders capable of inhibiting aggregation. They isol

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA