
Product Data Sheet

Product Name: Human CMV pp65 (495-503) trifluoroacetate salt

Cat. No.: GA24017

Chemical Properties

Cas. No. 153045-21-7

Formula $C_{42}H_{74}N_{10}O_{12}S$ M.Wt 943.18

Solubility 1mg/mL in water Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Immunodominant epitope derived from the human cytomegalovirus (HCMV) matrix protein pp65. It represents a useful tool for studying cell-mediated immunity against HCMV and might be used for the generation of HCMV-specific cytotoxic T-cells for adoptive

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA