
Product Data Sheet

Product Name: (d(CH²)¹,D-Ile²,Ile³,Arg⁴,Tyr⁵)-Vasopressin

Cat. No.: GA24028

Chemical Properties

Cas. No. 114438-90-3

Formula M.Wt 1178.47

Solubility Storage < -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

AVP antagonist specifically targeting the rat V2 receptor. Data obtained by Cotte et al. indicate that the type of amino acid in position 2 controls the rat versus human selectivity in AVP analogs.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA