
Product Data Sheet

Product Name: 2-Chlorotrityl chloride resin (100-200 mesh, 1.5-1.9 mmol/g)

Cat. No.: GA24057

Chemical Properties

Cas. No. 934816-82-7

Formula M.Wt

Solubility Storage < -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Example ADS

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA