
Product Data Sheet

Product Name: For-Met-Ala-Ser-OH

Cat. No.: GA24110

Chemical Properties

Cas. No. 17351-32-5

Formula $C_{12}H_{21}N_3O_6S$

M.Wt 335.38

Solubility DMSO : 100 mg/mL (298.17 mM; Need ultrasonic)

Storage < -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

fMAS was used as substrate in a continuous assay for peptide deformylase (PDF).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA