

---

**Product Data Sheet**

---

Product Name: H-Cys(Acm)-2-chlorotrityl resin (200-400 mesh, 0.4-0.7 mmol/g)

Cat. No.: GA24119

**Chemical Properties**

Cas. No.

Formula M.Wt

Solubility Storage < -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

**Background**

H-Cys(Acm)-2-CITrt resin is excellently suited for synthesizing peptides containing a C-terminal Cys(Acm).

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA