

---

## Product Data Sheet

---

Product Name: H-Cys(Trt)-2-chlorotrityl resin (200-400 mesh, 0.50-0.90 mmol/g)

Cat. No.: GA24121

### Chemical Properties

Cas. No.

Formula M.Wt

Solubility Storage < -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Peptides containing C-terminal cysteines can be obtained with practically no racemization by Fmoc-based SPPS using 2-chlorotrityl chloride resin.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA