

---

**Product Data Sheet**

---

Product Name: Histone H3 (21-44)

Cat. No.: GA24127

**Chemical Properties**

Cas. No. 1373516-71-2

Formula M.Wt 2505.91

Solubility Storage &lt; -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The histone fragment was used as substrate for the protein arginine methyltransferase PRMT4.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA