
Product Data Sheet

Product Name: HIV-1 tat Protein (47-57) amide

Cat. No.: GA24128

Chemical Properties

Cas. No. 697226-52-1

Formula M.Wt 1558.84

Solubility H₂O : ≥ 65 mg/mL (41.70 mM) Storage < -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

HIV-1 tat Protein (47-57) is a potent cell penetrating peptides.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA