
Product Data Sheet

Product Name: Nesfatin-1-Like Peptide (mouse)

Cat. No.: GA24149

Chemical Properties

Cas. No.

Formula M.Wt 8738.49

Solubility Storage < -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The insulinotropic peptide encoded by nucleobindin 1 upregulated preproinsulin mRNA expression and insulin secretion.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA