
Product Data Sheet

Product Name: RAD16-I
Cat. No.: GA24159

Chemical Properties

Cas. No. 289042-25-7

Formula M.Wt 1712.8

Solubility Storage < -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Self-assembling amphiphilic peptide used as scaffold in tissue engineering and 3-dimensional cell culture. Has also found use as carrier for drugs.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA