

---

**Product Data Sheet**

---

Product Name: Rink amide AM resin (200-400 mesh, 0.5-1.2 mmol/g)

Cat. No.: GA24160

**Chemical Properties**

Cas. No.

Formula M.Wt

Solubility Storage +5 ± 3 °C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-Rink amide-AM resin, a standard resin for the Fmoc-SPPS of peptide amides. Cleavage was performed with 50 % TFA in CH<sub>2</sub>Cl<sub>2</sub> or 95% aqueous TFA. Scavengers may be required.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA