
Product Data Sheet

Product Name: 5-Propynyl-dCTP

Cat. No.: GB20081

Chemical Properties

Cas. No.

Formula $C_{12}H_{18}N_3O_{13}P_3$ (free acid) M.Wt 505.2 (free acid)

Solubility Storage Store at -20°C or below

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA