
Product Data Sheet

Product Name: 2'-Azido-dATP

Cat. No.: GB20088

Chemical Properties

Cas. No.

Formula $C_{10}H_{15}N_8O_{12}P_3$ (free acid) M.Wt 532.2 (free acid)

Solubility Storage Store at -20°C or below

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

Sugar modified nucleoside triphosphates, increases target affinity and nuclease stability while reducing immune response.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA