

---

**Product Data Sheet**

---

Product Name: 5-Br-UTP  
Cat. No.: GB20102

**Chemical Properties**

Cas. No.

Formula  $C_9H_{14}N_2O_{15}P_3Br$  (free acid) M.Wt 563.0 (free acid)

Solubility Storage Store at -20°C or below

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

**Background**

Base modified ribonucleoside triphosphates, used in SELEX technique, provides thermal stability and nuclease resistance to aptamers.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA