
Product Data Sheet

Product Name: CP 615003 mesylate

Cat. No.: GB40156

Chemical Properties

Cas. No. 1259477-42-3

Formula $C_{20}H_{24}FN_3O_3 \cdot CH_4O_3S$

M.Wt 469.53

Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

A potent and subtype selective GABAA receptor partial agonist potentially useful in treating generalized anxiety disorder; the Mdr1 P-glycoprotein (P-gp) substrate

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA