
Product Data Sheet

Product Name: GSK 461364 analogue II

Cat. No.: GB40224

Chemical Properties

Cas. No. 929095-22-7

Formula C₂₇H₂₇F₃N₄O₃S

M.Wt

544.59

Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Potent Polo-like kinase (PLK) inhibitor, selective at PLK1 (IC₅₀: 2 nM) over PLK3 (IC₅₀: 270 nM)

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA