
Product Data Sheet

Product Name: LY 2955303 dihydrochloride

Cat. No.: GB40268

Chemical Properties

Cas. No. 2469014-43-3

Formula $C_{36}H_{42}N_4O_3 \cdot 2HCl$

M.Wt

651.67

Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Potent and selective retinoic acid receptor gamma (RAR γ) antagonist (Ki values >1700 nM, >2980 nM, and 1.1 nM for RAR α , RAR β , and RAR γ , respectively) for the treatment of osteoarthritis joint pain.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA