
Product Data Sheet

Product Name: TL 102 hydrobromide - DPAT, 5,6-Dihydroxy-

Cat. No.: GB40391

Chemical Properties

Cas. No. 62421-54-9

Formula $C_{16}H_{25}NO_2 \cdot HBr$

M.Wt

344.29

Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Dopamine receptor agonist

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA