
Product Data Sheet

Product Name: Dextran sulfate sodium salt (M.W 200000)

Cat. No.: GC10520

Chemical Properties

Cas. No. 9011-18-1

Formula $(C_6H_7Na_3O_{14}S_3)_n$ M.Wt Mr
~200,000

Solubility Water: 100 mg/mL clear, colorless to faintly greenish-yellow

Storage Store at 2-8°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA