
Product Data Sheet

Product Name: Supra Sieve GPG

Cat. No.: GC10986

Chemical Properties

Cas. No.

Formula M.Wt

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

analytical and preparative electrophoresis of DNA <1000bp PCR reactions, northern/southern blotting ligate and transform directly in melted gel no detectable inhibition of DNA ligase or the following restriction enzymes: BamH I, Hind II, EcoR I, Hind III, Pst I DNase, RNase, Protease: none detected store dry at room temperature, CAS 9012-36-6

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA