
Product Data Sheet

Product Name: DL-AP4 Sodium salt

Cat. No.: GC12960

Chemical Properties

Cas. No. 1263093-79-3

Chemical Name sodium (R)-hydrogen(3-amino-3-carboxypropyl)phosphonate

SMILES [O-][P@](O)(CC[C@H](C(O)=O)N)=O.[Na+]

Formula C₄H₉NNaO₅P M.Wt 205.08

Solubility <20.51mg/ml in Water Storage Desiccate at RT

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA