

---

**Product Data Sheet**

---

Product Name: Pridinol Methanesulfonate

Cat. No.: GC13661

**Chemical Properties**

Cas. No. 6856-31-1

Chemical Name 1,1-diphenyl-3-piperidin-1-ylpropan-1-ol;methanesulfonic acid

SMILES CS(=O)(=O)O.C1CCN(CC1)CCC(C2=CC=CC=C2)(C3=CC=CC=C3)OFormula  $C_{21}H_{29}NO_4S$  M.Wt 391.52Solubility  $\geq 16.95\text{mg/mL}$  in DMSO Storage Store at  $-20^{\circ}\text{C}$ General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}\text{C}$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}\text{C}$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Pridinol Methanesulfonate Salt is used as an anticholinergic and as an antiparkinsonian agent.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA